

Voor iedereen, die geen enkele uitdaging uit de weg gaat

INNOVATIEF TEAR-INFUSED DESIGN NU BESCHIKBAAR VOOR ASTIGMATISME

Met het grootste parameterbereik van alle
bestaande zachte torische daglenzen.^{1,2}



NEW
ACUVUE® 1-Day
oasys WITH HydraLuxe™
BRAND CONTACT LENSES
ASTIGMATISM


ACUVUE OASYS® 1-Day for ASTIGMATISM werkt samen met de natuurlijke traanfilm van het oog en biedt de hele dag lang uitzonderlijk comfort en consistent, duidelijk en stabiel zicht.^{1,3}

47%

VAN DE CONTACTLENSDRAGERS HEEFT EEN CILINDER VAN $\geq 0.75D$ IN MINSTENS ÉÉN OOG⁴


- Voor deze dragers kan elke hoofd- of oogbeweging de kwaliteit van het zicht beïnvloeden³
- Een juiste correctie kan de visuele kwaliteit van het leven van contactlensdragers enorm verbeteren⁵

Twee gecombineerde innovaties helpen bij dragers met astigmatisme het gevoel van vermoeide ogen gedurende de dag te voorkomen⁶



Tear-Infused Design

Een uniek design dat werkt door te integreren met de eigen traanfilm van de drager.¹



Eyelid stabilised design

Werkt op een natuurlijke manier met de oogleden om de lens in de juiste positie te houden.⁷



**UITZONDERLIJK
COMFORT EN
CONSISTENT,
DUIDELIJK EN
STABIEL ZICHT, DE
HELE DAG LANG^{1,3}**

De voordelen van ACUVUE OASYS® 1-Day for ASTIGMATISM

Voor uw praktijk...

99%

**SUCCES BIJ EERSTE
AANSPANNING**

Het grootste parameterbereik
van alle bestaande zachte torische daglenzen²

Voor uw klanten...

100%

ZICHTSCHERPTE⁶

9 op de 10 klanten ervaren een scherp en
duidelijke zicht, de hele dag lang⁶

Het innovatieve Tear-infused design is nu ook verkrijgbaar voor astigmatisme, met het grootste parameterbereik van alle bestaande zachte torische daglenzen.⁴

PARAMETERS										
Sferische sterktes	LAGE MIN (25) 0.00D t/m -6.00D (per 0.25D)				HOGE MIN (6) -6.50D t/m -9.00D (per 0.50D)			PLUS (16) +0.25D t/m +4.00D (per 0.25D)		
Assen/Cilinders	-0.75	-1.25	-1.75	-2.25	-0.75	-1.25	-1.75	-0.75	-1.25	-1.75
10°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
20°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
30°	✓	✓	✓							
40°	✓	✓	✓							
50°	✓	✓	✓							
60°	✓	✓	✓							
70°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
80°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
90°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
100°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
110°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
120°	✓	✓	✓							
130°	✓	✓	✓							
140°	✓	✓	✓							
150°	✓	✓	✓							
160°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
170°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
180°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓



NEW
ACUVUE® 1-Day
oasys WITH HydraLuxe™
 BRAND CONTACT LENSES
ASTIGMATISM

PRODUCTSPECIFICATIES

Materiaal	senofilcon A
Lenstechnologie	Polyvinyl pyrrolidone (PVP) met Hydraluxe™ technologie
Watergehalte	38%
Basiscurve	8.5 mm
Diameter	14.3 mm
Middendikte	0.080 mm (-3.00D)
Zuurstoftransmissie (Dk/t) (rand- en grensgecorrigeerd)*	129 x 10 ⁻⁹ (-3.00D)
Oxygen flux** (% beschikbaar voor centrale cornea)	98% (geopend oog)
UV-filter***	Klasse 1, 99% UVB en 90% UVA
Hanteerbaarheid	lichtblauwe handling tint
Vervangingsschema	Eenmalig gebruik: dagelijks vervangen
Verpakking	30 lenzen

STERKTEBEREIK

0.00D t/m -6.00D (per 0.25D)	Cilinders: -0.75D, -1.25D en -1.75D Assen: 10°-180° (per 10°)
0.00D t/m -6.00D (per 0.25D)	Cilinder: -2.25D Assen: 70°-110° en 160°- 20° (per 10°)
-6.50D t/m 9.00D (per 0.50D)	Cilinders: -0.75D, -1.25D, en -1.75D Assen: 70°-110° en 160°- 20° (per 10°)
+0.25D t/m +4.00D (per 0.25D)	Cilinders: -0.75D, -1.25D, en -1.75D Assen: 70°-110° en 160°- 20° (per 10°)

1. JJVC data on file 2016. Hydraluxe™ Technology definition. 2. JJVC data on file 2016. Proportion astigmats accommodated with ACUVUE® Brand Contact Lenses for ASTIGMATISM. 3. JJVC data on file 2017. Clinical and subjective performance of ACUVUE OASYS® 1-Day Brand Contact Lenses with Hydraluxe™ Technology for ASTIGMATISM. 4. Young G, Sulley A, Hunt C. Prevalence of astigmatism in relation to soft contact lens fitting. Eye Contact Lens. 2011;37(1):20-25. Of the 47.4% of patients, nearly half exhibited significant astigmatism in one eye only. 5. Nichols J et al. A comparison of toric and spherical soft contact lenses on visual quality of life and ease of fitting in astigmatic patients. BCLA/NCC conference 2016. In a clinical trial of habitual contact lens wearers with low levels of astigmatism as assessed using the NEI-RQL validated instrument, correction with 1-DAY ACUVUE® MOIST for ASTIGMATISM showed improved visual quality of life scores compared to correction with 1-DAY ACUVUE® MOIST after 1 week of wear. 6. Straker B, Hamada W, Sulley A, Olivares G. Fitting performance and efficiency with a low silicone hydrogel daily disposable toric contact lens. Poster presented at: GSLS Conference 2017. 100% of eyes experienced acceptable lens fit at fitting, with 99% experiencing acceptable fit at 1-week follow-up (321/324). Rx range: -1.50 to -4.00DS inclusive, -0.75 to 1.50DC inclusive and axes within 180 } 15° or 90 } 15°. Acceptable fit defined as absolute lens misrotation ≤20 degrees, and lens stability with blink ≤5 degrees. In a clinical trial, 97% of patients achieved a monocular VA 20/20 or better at the fitting visit with 100% achieving 20/25 or better. 7. JJVC data on file 2016. Design enhancements and resultant benefits of ACUVUE OASYS® Brand Contact Lenses 1-Day with Hydraluxe™ Technology for ASTIGMATISM.